
## *Conelius.* Quest Elite 4.5 L Pouch Holders & Shelf Conversion

Upon receipt of the retrofit kit containing 1 x new shelf and 2 x new 4.5 liter pouch holders, proceed on replacing the parts **ONLY after your current stock of old style juice concentrate pouches has depleted** and <u>preferably in combination with</u> the weekly cleaning. Time required to perform the conversion: 10 minutes ( excl. cabinet cleaning ) Tools required : None ( except for calibrating the Ratios.) Scan the QR-Code for viewing these instructions in other languages



4	Bag Loading Instructions - New Pouches: Load the new pouch with the large connector into the Large Slot. Make sure to install the connector as shown. Failure to install the connector in the right way will cause issues connecting the silicone tube to the mixing chamber, concentrate leaks or overall poor quality of the dispensed drinks.	
5	<b>Bag Loading Instructions – Old Pouches:</b> <b>!</b> You will note that the new pouch holders still have the smaller slot opposite the now required larger slot. This is in the event that the old style pouches may be used by way of a contingency back-up if ever required.	
6	Install the pouch holder onto the shelf and connect the silicone tube to the pump, close the pump door. Next, remove the red stoppers and put these aside in a dry and clean place. Install the mixing chambers and nozzles.	
7	Repeat the above for the second pouch holder / flavour. For details, please refer to the Quest Elite Operators Manual and PM cards, instruction videos or other.	
8	Check the Ratio. If out of specification, please follow your standard protocol for adjusting the ratio.	

٦

Г

1